

Waters® P200X PM Kit List

The following parts are included in the WatersP200X Preventive Maintenance Kit - 205001740. Keep this with the Waters *System Information Set*.

Part Number	Quantity	Item	Description/Location
289011368	2	ASSY, 200 CFM SAPPIRE / PISTON FRONT	Used when replacing the pistons
01033	2	Seal, P-200A, GFPM, HVY-Med Spring	Main seals in the pump
01032	2	Ring, P-200 Piston Seal Peek Backup	Backup rings for the seals
410006896	1	FILTER MEDIA, 80mm X 80mm X 5mm	Filter on the air Inlet of the pump
289012441	4	ASSY, P200X CHECK VALVE WITH O RINGS	Prevents backflow of fluid

