


## Symantec Endpoint Protection v12.1 Screenshots

System: Xevo TQ-S , Windows 7, software test lab, Waters Wilmslow, 1st August 2016



Change Settings - Symantec Endpoint Protection

## Change Settings

Help

Status

- Scan for Threats
- Change Settings**
- View Quarantine
- View Logs
- LiveUpdate...

You can change the following settings:

- 
**Virus and Spyware Protection**  
 Protects against viruses, malware, and spyware  
[Configure Settings](#)
- 
**Proactive Threat Protection**  
 Provides zero-day protection against unknown threats  
[Configure Settings](#)
- 
**Network Threat Protection**  
 Protects against Web and network threats  
[Configure Settings](#)
- 
**Exceptions**  
 Exceptions provide the capability to exclude files, folders, and processes from being scanned  
[Configure Settings](#)
- 
**Client Management**  
 Provides functionality to manage this client  
[Configure Settings](#)




Scan for threats - Symantec Endpoint Protection

## Scan for threats


Help

Status

- Scan for Threats**
- Change Settings
- View Quarantine
- View Logs
- LiveUpdate...



**Active Scan**  
Scan only the most commonly infected areas.  
Duration: **1-2 minutes**  
[Run Active Scan](#)




**Full Scan**  
Scans your entire computer and can run in the background.  
Duration: **30-120 minutes**  
[Run Full Scan](#)


**Scans**  
Scans configured for this computer. [View Global Scan Settings](#)

[Create a New Scan](#)

Scan Name	Enabled	Type	When to Scan	Last Scan
Monthly Scheduled Scan	Yes	Full Scan	Monthly	July 01, 2016 10:37 AM
Weekly Scheduled Scan	Yes	Active Scan	Weekly	July 29, 2016 1:01 PM




## 1. Virus and Spyware Protection settings




Click View List: no exceptions defined




click Notifications:



Click Actions:



click Advanced:



Virus and Spyware Protection Settings

- Outlook Auto-Protect
- Notes Auto-Protect
- Early Launch Anti-Malware
- Global Settings
- Auto-Protect
- Download Insight
- Internet Email Auto-Protect

Enable Download Insight to detect potential risks based on file reputation

[What is File Reputation?](#)

Actions

Notifications

Specify download sensitivity level

Select level:

- 9 (Maximum)
- 8
- 7
- 6 (High)
- 5 (Typical)
- 4
- 3
- 2
- 1 (Minimum)

Level 5 (Typical)

Allows only files that have a good or unknown reputation. Some files are considered malicious, and some files are considered unproven. The number of false positive detections is low.

Also detect files as malicious based on their use in the Symantec Community

Files with: 5 or fewer users

Files known by users for: 2 or fewer days

Intranet Settings


Automatically trust any file downloaded from an intranet website

OK


Cancel

Help


click Actions:









click Actions;




click Notifications




click Advanced:








click Actions:




Click Notifications:




click on Advanced:





Click on Actions



Click on Notifications;






click on Advanced:









## 2. Proactive Threat Protection settings






### 3. Network Threat Protection settings







#### 4. Exceptions



#### 5. Client Management Settings



Proxy Options:



The image shows a 'Proxy Server Settings' dialog box with a close button (X) in the top right corner. It features a tab labeled 'HTTP or HTTPS'. The main content area is titled 'Proxy Configuration' and contains the following elements:

- A text instruction: "Select the proxy server option to use for all external communication except LiveUpdate."
- Three radio button options:
  - I do not use a proxy server
  - I want to use the proxy server specified by my client browser (default)
  - I want to use custom proxy settings
- Input fields for custom settings:
  - Host proxy: [text box]
  - Port: [text box]
  - Authentication required
  - User name: [text box]
  - Password: [text box]



At the bottom of the dialog are three buttons: 'OK', 'Cancel', and 'Help'.

## Reboot Options:

The image shows a Windows dialog box titled "Reboot Options". It is divided into four sections: "Restart Computer", "Restart Type", "Restart Message", and "Other Options".

- Restart Computer:** Includes radio buttons for "Immediately" and "At this time:". The "At this time:" option is selected, with a dropdown menu showing "Today" and a time field set to "6:00:00 AM". There is also a checked checkbox for "Randomize start time" and a "Randomization Hours" spinner set to "3".
- Restart Type:** Includes radio buttons for "Do not restart (not recommended)", "No prompt", "Prompt with a countdown of:", and "Prompt and allow snooze up to:". The "Prompt and allow snooze up to:" option is selected, with a spinner set to "5" and the text "minutes" and "Scheduled Time".
- Restart Message:** A text box containing the message: "A security risk remediation or a new content download requires this computer to restart."
- Other Options:** Includes checked checkboxes for "Hard restart" (with a warning: "Third party applications will not be able to prevent this restart. Unsaved user data could be lost.") and "Restart immediately if the user is not logged in" (with a warning: "If the user is not logged in when the restart request is sent, this option forces an immediate restart and overrides other pending restart actions.").

At the bottom of the dialog are three buttons: "OK", "Cancel", and "Help".





click Configure Proxy Options:

LiveUpdate Proxy Server Settings

HTTP or HTTPS FTP

HTTP or HTTPS Proxy Configuration

Select the type of configuration for communications

I do not use a proxy server for HTTP or HTTPS

I want to use my Windows Internet Options proxy settings (default)

I want to customize my HTTP or HTTPS settings

Host proxy:

HTTP Port:  HTTPS Port:

Authentication required

Authentication Type

Basic Authentication

User name:

Password:

NT LAN Manager Authentication

User name:

Password:

Domain name:

OK Cancel Help

LiveUpdate Proxy Server Settings



HTTP or HTTPS

FTP

FTP Proxy configuration

Select the type of configuration for communications

- I do not use a proxy server for FTP
- I want to use my Windows Internet Options FTP settings (default)
- I want to customize my FTP settings

FTP host proxy:

Port:

OK

Cancel

Help

Client Management Settings


General | Tamper Protection | LiveUpdate | Submissions

Submissions help Symantec address threats in the fight against digital cybercrime.


Let this computer automatically forward selected anonymous security information to Symantec

Select additional types of information to submit:

- File reputation
- Antivirus detections
- Antivirus advanced heuristic detections
- SONAR detections
- SONAR heuristic

 If you disable any submission types, you prevent Symantec Endpoint Protection from providing the best protection for your client computers.

Allow Insight lookups for threat detection (recommended)

 If you disable this option, you prevent Download Insight from detecting threats on your client computers. You also limit the protection that is provided by virus and spyware scans and SONAR. The rate of false positives might increase.

OK

Cancel

Help