

Empower QuickStart and Empower QuickStart Network

Chromatography data software for your workstation or small network

HELP YOUR LABORATORY OPERATE MORE EFFICIENTLY

Affordable and easy-to-use, Empower™ QuickStart and Empower QuickStart Network offer robust features that strike a balance of simplicity, performance, and reliability. Providing the perfect software and system solution for any cost-conscious organization with a limited budget and limited chromatography experience. Available as workstation or small network, Empower QuickStart enables your laboratory to operate more efficiently.

INCREASE PRODUCTIVITY AND DECREASE TRAINING COSTS

Get started quickly with a user-friendly interface that requires minimal training, and is easy to set up, learn, and operate, allowing you to acquire, process, and report LC analyses — converting high-quality and accurate results into the answers you require. It simplifies the complete management of your chromatographic process and data, providing reliable information so you can make informed and accurate decisions quickly — all the while offering a high level of success with low cost of ownership.




Figure 1. The Empower QuickStart Interface uses a navigation bar to guide you through the software and all tasks are carried out in one window.

EASILY ACCESSIBLE

Use a single window to access simplified tools for:

- Chromatographic instrument control
- Chromatographic data acquisition
- Interactive data processing and data review
- Comprehensive reporting
- Management of data and chromatography system

ENSURE DATA TRACEABILITY AND DATA INTEGRITY

Having confidence in the quality of your data and in your results is critical. Automatically manage data through several levels of safeguards to ensure trustworthy and reliable records. Empower QuickStart makes it easy to review and approve results, search for data, and verify records for authenticity. It simplifies the complete management of your chromatographic workflow and data, providing reliable information so you can make informed and accurate decisions efficiently. When your laboratory size or feature requirements grow, all data collected can be upgraded to the global standard Empower Software.

HELP WHEN YOU NEED IT

Setting up new acquisitions and processing chromatographic results can sometimes be challenging – especially when the task is unfamiliar. Complete tasks quickly and more accurately with Show Me... Help, a guided workflow wizard.




Figure 2. Show Me... Help provides interactive help in carrying out common tasks.

FLEXIBLE OPTIONS TO FIT YOUR NEEDS

Choose either Empower QuickStart Workstation or Empower QuickStart Network – both system solutions are flexible and built to suit your needs:

	Empower QuickStart	Empower QuickStart Network
Licenses	5	5 (up to 20)
Systems	2	Connects up to 10
Interfaces	QuickStart	QuickStart
Instrument control	Waters Alliance™ 2489, and 2998 Detectors	Multi-vendor control
e-SAT/IN	Optional	Optional
Advanced options:		
System suitability	Included	Included
GPC	Optional	Optional

WATERS GLOBAL SERVICES

Delivering world renowned services and support

Waters Service and Support offerings are tailored to optimize your laboratory productivity while addressing your budget realities. Our offerings help maintain system peak performance, minimize down time, address scientific application challenges, and support stringent compliance requirements. As your services and support provider, we are committed to the success of your laboratory and business.

Waters quality support and consultative services ensure your success wherever your laboratory is located in the world.



Waters has consecutively earned the ACE Award since 2001 for providing best-in-class technical knowledge, issue resolution, and process support.

Waters offers comprehensive and cost-effective Software Support/ Maintenance and Qualification Plans to maximize the value of your investment. To learn more, please contact your local Waters office.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Waters, The Science of What's Possible, Empower, and Alliance are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.