



Full Install

FR: Feature Release
 SR: Service Release
 HF: Hot Fix

Beginning with Empower 3.6.0
 the naming convention for
 Empower release has changed to
 XX.YY.ZZ.QQ where
 XX= Major Release Version
 YY= Feature Release Version
 ZZ= Service Release Version
 QQ= Hotfix Release Version

Empower 3 Base
 Oracle 11.2.0.1.0
 DB 7.00.00.00

Empower 3 SR1
 Oracle 11.2.0.1.0
 DB 7.00.00.99

Empower 3 FR1
 Oracle 11.2.0.2 patch set 11
 DB 7.10.00.00

Empower 3 FR2
 Empower 3 FR2 HF1
 Oracle 11.2.0.2 Win
 Oracle 11.2.0.3.6 RedHat & Solaris
 DB 7.20.00.00

Empower 3 SR2
 Empower 3 SR2 HF1
 Empower 3 SR2 HF2
 Empower 3 SR2 HF3
 Empower 3 SR2 HF4
 Oracle 11.2.0.2 Win
 Oracle 11.2.0.3.6 RedHat & Solaris
 DB 7.21.00.00

Empower 3 FR3
 Empower 3 FR3 HF1
 Oracle 11.2.0.4
 DB 7.30.00.00

Empower 3 FR4
 Oracle 12.1.0.2
 DB 7.40.00.00

Empower 3 SR3
 Empower 3 SR3 HF1
 Empower 3 SR3 HF2
 Empower 3 SR3 HF3
 Oracle 12.1.0.2
 DB 7.41.00.00

Empower 3 FR5
 Empower 3 FR5 HF1
 Empower 3 SR4
 Oracle 18.3.0.0 - 32Bit Client
 Oracle 18.3.0.0 - 64Bit DB RedHat
 Oracle 18.4.0.0 - 64Bit DB Windows
 DB 7.50.00.00

Empower 3 SR5
 Empower 3 SR5 HF1
 Oracle 18.3.0.0 - 32Bit Client
 Oracle 18.3.0.0 - 64Bit DB RedHat
 Oracle 18.4.0.0 - 64Bit DB Windows
 DB 7.51.00.00

Empower 3.6.0
 Oracle 19.3.0.0 - 32Bit Client
 Oracle 19.7.0.0 - 64Bit DB RedHat
 Oracle 19.7.0.0 - 64Bit DB Windows
 DB 7.60.00.00

Empower 3.6.1
 Oracle 19.3.0.0 - 32Bit Client
 Oracle 19.7.0.0 - 64Bit DB RedHat
 Oracle 19.7.0.0 - 64Bit DB Windows
 DB 7.61.00.00

Empower 3.7.0
 Oracle 19.3.0.0 - 32Bit Client
 Oracle 19.7.0.0 - 64Bit DB RedHat
 Oracle 19.7.0.0 - 64Bit DB Windows
 DB 7.70.00.00

Empower 3.8.0
 Oracle 19.3.0.0 - 32Bit Client
 Oracle 19.14.0.0 - 64Bit DB RedHat
 Oracle 19.14.0.0 - 64Bit DB Windows
 DB 7.80.00.00

